

# PARTIAL HAND

## CLINICIAN STEPS

### 1. PATIENT MEASUREMENT

#### Cast patient

- Hand must be in proper position – fingers flat, wrist neutral. Mold must extend from the distal tip of fingers to a minimum of 2" past the styloid (up to the cubital fold is best).
- Use only Body Double or P386 (Gray) impression material
- Refer to resource materials at Partial Hand Solutions for proper position and casting technique.

#### Complete the Silicone Interfacing order form

- Refer to Ottobock color chart for skin-tone silicone color options. Black and solid primary colors are also available.
- Though not required, photos, measurements, and tracings of patient are encouraged.

#### Order parts (M-Fingers and TITAN fingers)

- Fingers are separate components
- Vendor options: Hanger Direct, College Park, LTI, or the Partial Hand Solutions website.
- Componentry costs are separate from fabrication costs.
- Ship cast, components and order form to:  
Hanger Fabrication Network / Silicone Interfacing Department  
10 County Line Drive, Unit E  
Cromwell, CT 06416

### CHECK SOCKET FABRICATION

On receipt on cast components and order form:

- SI Department will contact your office to review details and provide a timeline for fitting.
- Standard production times are approximately 2-6 weeks for a check fitting, depending on scheduling and complexity of the device.
- On check socket fabrication completion, we will ship to your office.

### 2. CHECK SOCKET FITTING

- Fittings are required for most devices – please fit and align parts to patient specifications.
- 'Straight to finish' devices are not recommended, the office will be held financially responsible for any additional adjustments or remakes.
- Refer to resource materials at Partial Hand Solutions for proper alignment of M-Fingers and Titan fingers.
- Return check socket and components to:  
Hanger Fabrication Network / Silicone Interfacing Department  
10 County Line Drive, Unit E  
Cromwell, CT 06416

### 3. DEFINITIVE DEVICE FABRICATION

- We will contact your office with a definitive ship date.
- Standard production times are approximately 2-4 weeks for completing a definitive device from fitting, depending on production schedule and device complexity.
- On completion, we will ship the definitive device to your office for patient delivery.

