



MEERKAT²

PEDIATRIC LOW TONE SMO

INDICATIONS

- Hypotonia
- Flexible Pes planus, moderate to severe
- Ankle instability, moderate to severe
- No cast required for children < 80 lbs

PATIENT-CENTERED

- Tailored to each patients needs through data-controlled designs
- Timely delivery (5 days)
- Kid-friendly packaging, AFO socks and stuffed animal
- Guaranteed fit*

UNMATCHED EXPERIENCE

- Developed by clinicians, engineers and technicians
- Regulated by proprietary big data algorithm
- Directed by National Standards of Care and Hypotonia Clinical Practice Guidelines

HOW TO ORDER



Download and complete order form: www.hangerfabrication.com



Email form to: HFN_Orlando@hanger.com



5 day turn-around time

CONTRAINDICATIONS

- Clubfoot, forefoot deformities, high tone

FEATURES

- Sustentaculum tali modification and internal medial heel post

* GUARANTEED FIT: Simply provide us with a *complete* order form, we guarantee the fit.

Note critical measurements that drive fit are: 1] length from heel to 1st met apex and 2] ML at met heads.